





Sound device setup

Before start of Clinic or before creating new documents on the G2 system, it is **very** important that Authors ensure that's the Speech Mike is set as the default recording device in G2.

- Please click on the cog wheel located towards the top right corner of the G2 home screen to open the settings menu.



- Go to the recording tab in the G2 menu and ensure that the speech Mike is listed as the recording device as in the illustration below.



- Perform microphone calibration (if necessary) by clicking Start Microphone Calibration and following the on screen instructions.
- When finished, click “Save & Close” in the top left-hand corner.

