



















Printing and Reports	
Select the Home tab	
With a Patient chart open click Print	
You can then tick which document(s) you wish to print from the selection	
With a Patient open you will be able to print documents within their ICCA record	
Click Preview to open a preview of the document prior to printing	
Change any printer settings at the top of the window	
Click OK	

Printing and Reports	
Some reports will have extra options that you can select at this stage	 
Review these settings for each document by selecting each one on the left-hand side by clicking Next	 
Click Print	

Handover Documents	
Select the Home tab	
With a Unit	
You can then tick the census report from the selection	
With the Census screen of a unit selected you can print a Census handover document including patient allergies and alerts	
Click Preview to open a preview of the document prior to printing	
Change any printer settings at the top of the window	
Click OK	
Alternatively:	
You can print a more detailed handover document via the Xlinks tab	
Click the Reports Page button	