












Saving Data and Data Indicators/Icons	
<p>Save the Patient record using the <b>Save Chart</b> button in the top-left of the chart tab</p> <p>OR</p> <p>Use <b>CTRL+S</b> on the keyboard</p>	
<p><b>Important</b> - Save regularly, if you do not save or interact with ICCA within 9 minutes the system will log you out and changes will be lost</p>	
<p>Black text</p> 	<p>Saved Data</p>
<p>Pink Text</p> 	<p>Validated (input data) but unsaved – Click <b>Save Chart</b> to turn this black</p> 
<p>Green Text</p> 	<p>Unvalidated (previously entered data in the same field – needs accepting) – Use <b>Accept Defaults</b> to accept either just the one selection or the whole form/chart that you're viewing</p> 
<p>Red Text</p>	<p>Conflicting data (saved or unsaved) i.e. two people trying to input data into the same cell/form</p>
<p>Strikethrough Data</p> 	<p>Deleted or Discontinued Row (Unsaved) – Click <b>Save Chart</b> to save the delete or discontinuation</p> 
<p>Greyed out text</p> 	<p>A saved discontinued row or unavailable/disabled feature or function</p>
<p>An Ellipsis “...”</p>	<p>Data has been truncated to fit the cell/field</p>

Saving Data and Data Indicators/Icons	
<p>Black triangle in the lower-right corner</p> 	<p>Data has been changed</p>
<p>Red triangle in the top-right</p> 	<p>Remarks have been added in the properties – For lab results this also means that the lab result does not fit the cell and you must open the remarks to view</p>